

essilor

Product Range



September 2025

Table of Contents



Contents

How To Use This Guide

How To Use This Guide	3
-----------------------------	---

Progressive Lenses

Varilux® XR series™	4
Varilux® Physio® Extensee™	6
Varilux® Comfort Max	9

Extended Focus Lenses

Varilux® Digitime™	12
--------------------------	----

Essilor Essentials

Essilor Essentials™ Assent	13
Essilor Essentials™ Accord	15

Digital Eyestrain Lenses

Eyezen™ Boost	17
---------------------	----

New Generation of SV Lenses

Eyezen™ Start	18
---------------------	----

SV 360 Fit

Aspheric	21
----------------	----

Single Vision Lenses

Aspheric	22
Spherical	23

Sun Lenses

Solarflair Progressive	25
Solarflair SV	26

Myopia Control Lenses

Stellest™	27
Myopilux®	28

Wrap Lenses

Varilux® Stylistic™	29
Single Vision Stylistic™	30

Extended Focus Lenses

Interview®	31
------------------	----

Bifocal & Trifocal Lenses

Bifocal & Trifocal	32
--------------------------	----

Special Lenses Range

Exceptio Special Lenses	33
-------------------------------	----

Stock Lenses

Stock Index 1.74	34
------------------------	----

Stock Lenses

Stock Index 1.67	35
Stock Index 1.60	36
Stock Index 1.59	38
Stock Index 1.50	40

Varilux XR

Base Chart 1.50	43
Base Chart 1.59	44
Base Chart 1.60	45
Base Chart 1.67 and 1.74	46

Lens Engravings


Lens Engravings	47
-----------------------	----

Lens Technology

Crizal® Range	48
Eye Protect System™	49
Blue UV Capture™	50
Transitions® GEN S™	51
Transitions® XTRActive® New Generation	52
Transitions® XTRActive® Polarized	52
Transitions® Style Mirrors	53

Sunwear

Essilor® Sun Solutions	54
Xperio® Tints	55
Essilor® Sun Solutions Availability	56



Legend

Icon	Product
Clear	Clear / Standard Lens
EyeProtect System™	Eye Protect System™
BlueUV Capture™	Blue UV Capture™
Transitions® Gen S™	Transitions® GEN S™ (Availability; see colour icons)
	Grey
	Brown
	Graphite Green
	Sapphire
	Amethyst
	Amber
	Emerald
	Ruby
XTRACTIVE® NEW GENERATION	Transitions™ XTRActive® New Generation Lens (Available in Grey, Brown & Graphite Green)
XTRACTIVE® POLARIZED™	Transitions™ XTRActive® Polarized Lens (Available in Grey)
Transitions® Drivewear®	Transitions™ Drivewear® (Available in Copper only)
Xperio®	Xperio® Polarised Lens (Availability; see colour icons)

Varilux® XR series™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire™^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

*Note: Varilux® XR series™ is only available with Crizal® and Optifog® coatings.

Versions

	MFH ⁽³⁾	Parameters
Varilux® XR Pro	14-17mm	Personalised / Visiooffice only: NVB (App Only), Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance, Dominant Eye, Head cape & Head / Eye coefficient
Varilux® XR 4D	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights & Eye Rotation Centre distance, Frame Wrap, Pantoscopic Tilt, Dominant Eye and Head Cape Angle
Varilux® XR Track	14-17mm	Personalised: NVB (App Only), Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® XR Fit	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® XR Short	14mm	Short corridor: Addition up to 3.50 D. Monocular Pupillary Distances and Fitting Heights
Varilux® XR Series™	17mm	Monocular Pupillary Distances and Fitting Heights

Range

1.74	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions	
Clear Transitions Gen[®]S EyeProtect System...	Varilux [®] XR series™						
	60	-14.00	+13.00	-6.00	0.75 to 4.00	Combined 6.50 DΔ	Regular Short Fit Track 4D Pro
	65	-13.00	+13.00	-6.00			
	70	-12.00	+10.00	-6.00			
	75	-11.00	+8.00	-5.00			
80	-10.50	+6.00	-4.00				

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions	
Clear Transitions Gen[®]S XTRACTIVE™ NEW GENERATION EyeProtect System...	Varilux [®] XR series™						
	60	-12.00	+9.00	-6.00	0.75 to 4.00	Combined 6.50 DΔ	Regular Short Fit Track 4D Pro
	65	-12.00	+9.00	-6.00			
	70	-11.00	+9.00	-6.00			
	75	-10.00	+8.00	-5.00			
80	-9.00	+6.00	-4.00				

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions	
Clear Transitions Gen[®]S XTRACTIVE™ NEW GENERATION EyeProtect System..	Varilux® XR series™						
	60	-12.00	+8.00	-6.00	0.75 to 4.00	Combined 6.00 DΔ	Regular Short Fit Track 4D Pro
	65	-12.00	+8.00	-6.00			
	70	-10.00	+8.00	-6.00			
	75	-9.00	+8.00	-5.00			
	80	-8.00	+6.00	-4.00			
	60	-12.00	+6.00	-6.00			
	65	-12.00	+6.00	-6.00			
	70	-10.00	+6.00	-6.00			
	75	-9.00	+6.00	-5.00			
80	-8.00	+6.00	-4.00				
85	-7.00	+4.00	-3.00				
Xperio.®							

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating not available in 1.67 and 1.74 indexes.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Varilux® XR series™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Prevencia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

*Note: Varilux® XR series™ is only available with Crizal® and Optifog® coatings.

Versions

	MFH ⁽³⁾	Parameters
Varilux® XR Pro	14-17mm	Personalised / Visioffice only: NVB (App Only), Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance, Dominant Eye, Head cape & Head / Eye coefficient
Varilux® XR 4D	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights & Eye Rotation Centre distance, Frame Wrap, Pantoscopic Tilt, Dominant Eye and Head Cape Angle
Varilux® XR Track	14-17mm	Personalised: NVB (App Only), Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® XR Fit	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® XR Short	14mm	Short corridor: Addition up to 3.50 D. Monocular Pupillary Distances and Fitting Heights
Varilux® XR Series™	17mm	Monocular Pupillary Distances and Fitting Heights

Range

1.59*	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Varilux® XR Series™						
Clear	60	-10.00	+6.00	-6.00	0.75 to 4.00	Combined 6.00 DΔ
Transitions Gen®S	65	-10.00	+6.00	-6.00		
EyeProtect System™	70	-9.50	+6.00	-6.00		
Xperio.	75	-8.50	+6.00	-5.00		
	80	-7.50	+5.00	-4.00		
						Regular Short Fit Track 4D Pro

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁵⁾	Versions
Varilux® XR Series™						
Clear	60	-10.00	+6.00	-6.00	0.75 to 4.00	Combined 5.00 DΔ
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-7.50	+6.00	-5.00		
	80	-6.00	+5.00	-4.00		
	85	-4.00	+3.50	-3.00		
Transitions Gen®S	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-7.50	+6.00	-5.00		
	80	-6.00	+5.00	-4.00		
XTRACTIVE® NEW GENERATION	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-7.50	+6.00	-5.00		
	80	-6.00	+5.00	-4.00		
EyeProtect System™	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-8.50	+6.00	-6.00		
	75	-7.00	+6.00	-5.00		
	80	-6.00	+5.00	-4.00		
Xperio.	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-8.50	+6.00	-6.00		
	75	-7.00	+6.00	-5.00		
	80	-6.00	+5.00	-4.00		
						Regular Short Fit Track 4D Pro

(1) Crizal® Easy Pro coating not available in 1.67 and 1.74 indexes.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Varilux®

Physio® extensee™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire™^{HR}

Crizal®
Prevencia®

Crizal®⁽²⁾
Sun XProtect

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Physio® Extensee™	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Fit	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)
Varilux® Physio® Extensee™ Track	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Track Eyecode	14-17mm	Personalised / Visiooffice only: NVB (App Only), Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance

Range

1.74	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear EyeProtect System™	60	-18.00	+9.00	-6.00	0.75 to 4.00	Combined 6.50 DΔ Regular Short Fit Eyecode
	65	-15.00	+9.00	-6.00		
	70	-13.50	+8.00	-5.00		
	75	-10.00	+7.00	-4.00		
	80	-9.00	+6.00	-3.00		
Transitions Gen®S	60	-19.00	+13.00	-6.00	0.75 to 4.00	Combined 6.50 DΔ Regular Short Fit Eyecode
	65	-15.00	+13.00	-6.00		
	70	-12.00	+9.00	-6.00		
	75	-10.00	+8.00	-4.00		
	80	-9.00	+6.00	-3.00		

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	55 / 60	-14.00	+9.00	-6.00	1.00 to 3.00	Combined 6.50 DΔ Regular Short Fit Eyecode
	65	-12.00	+9.00	-6.00		
	70	-11.00	+9.00	-6.00		
	75	-10.00	+8.00	-5.00		
	80	-9.00	+6.00	-4.00		
Transitions Gen®S EyeProtect System™	60	-12.00	+9.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Regular Short Fit Eyecode
	65	-12.00	+9.00	-6.00		
	70	-11.00	+9.00	-6.00		
	75	-10.00	+8.00	-5.00		
	80	-9.00	+6.00	-4.00		
XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED	60	-14.00	+9.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Regular Short Fit Eyecode
	65	-12.00	+9.00	-6.00		
	70	-11.00	+9.00	-6.00		
	75	-10.00	+8.00	-5.00		
	80	-9.00	+6.00	-4.00		
Xperio®	60	-12.00	+9.00	-4.00	0.75 to 3.50	Combined 6.00 DΔ Regular Short Fit Eyecode
	65	-12.00	+9.00	-4.00		
	70	-11.00	+9.00	-4.00		
	75	-10.00	+8.00	-4.00		
	80	-9.00	+6.00	-4.00		

(1) Crizal® Easy Pro coating is not available in 1.67 and 1.74 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Available in Transitions® Grey only.

(5) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details.

Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter.

Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

Varilux®

Physio® extensee™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Prevencia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Physio® Extensee™	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Fit	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)
Varilux® Physio® Extensee™ Track	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Track Eyecode	14-17mm	Personalised / Visiooffice only: NVB (App Only), Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance

Range

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Varilux® Physio® Extensee™						
Clear Transitions Gen [®] S™	60 / 65	-12.00	+8.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ
	65 / 70	-10.00	+8.00	-6.00		
	70 / 75	-9.00	+8.00	-5.00		
	75 / 80	-8.00	+5.00	-4.00		
EyeProtect System™	60	-12.00	+8.00	-6.00		
	65	-12.00	+8.00	-6.00		
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
XTRACTIVE™ NEW GENERATION	60	-12.00	+8.00	-6.00		
	65	-12.00	+8.00	-6.00		
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
XTRACTIVE™ POLARIZED™	60	-12.00	+8.00	-6.00		
	65	-12.00	+8.00	-6.00		
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
Xperio.	60	-12.00	+8.00	-4.00		
	65	-12.00	+8.00	-4.00		
	70	-10.00	+8.00	-4.00		
	75	-9.00	+8.00	-4.00		
Xperio.	60	-12.00	+8.00	-4.00		
	65	-12.00	+8.00	-4.00		
	70	-10.00	+8.00	-4.00		
	75	-9.00	+8.00	-4.00		

1.59*	Ø	Range	Cyl Max	ADD	Prism ⁽⁵⁾	Versions
Varilux® Physio® Extensee™						
Clear Transitions Gen [®] S™	60 / 65	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ
	65 / 70	-9.50	+6.00	-6.00		
	70 / 75	-8.50	+6.00	-5.00		
	75 / 80	-7.50	+5.00	-4.00		
EyeProtect System™ XTRACTIVE™ POLARIZED™	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-9.50	+6.00	-6.00		
	75	-8.50	+6.00	-5.00		
Xperio.	60 / 65	-10.00	+6.00	-6.00		
	65 / 70	-9.50	+6.00	-5.00		
	70 / 75	-8.50	+3.50	-4.00		
	75 / 80	-7.50	+5.00	-4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Available in Transitions® Grey only.

(5) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details.

Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter.

Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.



Varilux®

Physio® extensee™

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Physio® Extensee™	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Physio® Extensee™ Fit	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)
Varilux® Physio® Extensee™ Track	14-17mm	Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Physio® Extensee™ Track Eyecode	14-17mm	Personalised / Visiooffice only: NVB (App Only), Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance

Range

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁵⁾	Versions
Clear 	60	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 5.00 DΔ
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-6.00	+6.00	-5.00		
	80	-4.00	+4.00	-4.00		
XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED®	60	-10.00	+6.00	-6.00		
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-8.00	+6.00	-5.00		
EyeProtect System™	60	-10.00	+7.00	-6.00		
	65	-10.00	+7.00	-6.00		
	70	-9.00	+7.00	-6.00		
	75	-7.50	+6.00	-5.00		
Xperio.	60	-10.00	+6.00	-4.00		
	65	-10.00	+6.00	-4.00		
	70	-8.00	+6.00	-4.00		
	75	-7.00	+4.00	-4.00		

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Available in Transitions® Grey only.

(5) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details. Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter. Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.

Varilux® Comfort Max



Coatings

Crizal®
Easy^{Pro}

Crizal®
RockTM

Crizal®
Sapphire^{HR}

Crizal®
Previncia[®]

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Comfort Max	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Fit	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Comfort Max Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.74	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear Transitions Gen®S EyeProtect System..	Varilux® Comfort Max					
	55	-20.00	+13.00	-6.00	0.75 to 3.50	Combined 6.50 DΔ
	60	-15.00	+13.00	-6.00		
	65	-12.00	+9.00	-6.00		
	70	-10.00	+8.00	-5.00		
	75	-9.00	+6.00	-5.00		
						Regular Short Fit Eyecode

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear Transitions Gen®S BlueUV Capture.. EyeProtect System..	Varilux® Comfort Max					
	55	-14.00	+9.00	-6.00	0.75 to 3.50	Combined 6.50 DΔ
	60	-12.00	+9.00	-6.00		
	65	-11.00	+9.00	-6.00		
	70	-10.00	+8.00	-5.00		
	75	-9.00	+6.00	-4.00		
XTRACTIVE®* NEW GENERATION XTRACTIVE® POLARIZED®	60	-14.00	+9.00	-6.00		
	65	-12.00	+9.00	-6.00		
	70	-11.00	+9.00	-6.00		
	75	-10.00	+8.00	-5.00		
	80	-9.00	+6.00	-4.00		
Xperio®	55	-12.00	+9.00	-4.00		
	60	-12.00	+9.00	-4.00		
	65	-11.00	+9.00	-4.00		
	70	-10.00	+8.00	-4.00		
	75	-9.00	+6.00	-4.00		
						Regular Short Fit Eyecode

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.74 and 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Varilux® Comfort Max



Coatings

**Crizal®
Easy^{Pro}**

**Crizal®
Rock™**

**Crizal®
Sapphire^{HR}**

**Crizal®
Prevencia®**

**Crizal®
Sun XProtect⁽²⁾**

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Comfort Max	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Fit	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Comfort Max Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear Transitions gen®S BlueUV Capture.. EyeProtect System..	Varilux® Comfort Max					
	55	-12.00	+6.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ
	60	-12.00	+8.00	-6.00		
	65	-10.00	+8.00	-6.00		
	70	-9.00	+8.00	-5.00		
	75	-8.00	+6.00	-4.00		
XTRACTIVE®* NEW GENERATION XTRACTIVE® POLARIZED*	60	-12.00	+8.00	-6.00		
	65	-12.00	+8.00	-6.00		
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
	80	-8.00	+6.00	-4.00		
Xperio.	55	-12.00	+8.00	-4.00		Regular Short Fit Eyecode
	60	-12.00	+8.00	-4.00		
	65	-10.00	+8.00	-4.00		
	70	-9.00	+8.00	-4.00		
	75	-8.00	+7.00	-4.00		

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear Transitions gen®S EyeProtect System..	Varilux® Comfort Max					
	55	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ
	60	-10.00	+6.00	-6.00		
	65	-9.50	+6.00	-6.00		
	70	-8.50	+6.00	-5.00		
	75	-7.50	+5.00	-4.00		
BlueUV Capture.. Xperio.	55	-10.00	+6.00	-6.00		
	60	-10.00	+6.00	-6.00		
	65	-9.50	+6.00	-6.00		
	70	-8.50	+6.00	-5.00		
	75	-7.25	+5.00	-4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Varilux® Comfort Max

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire™^{HR}

Crizal®
Prevencia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Varilux® Comfort Max	17mm	Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights
Varilux® Comfort Max Fit	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Varilux® Comfort Max Eyecode	14-17mm	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Varilux® Comfort Max					
	55	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 5.00 DΔ
	60	-10.00	+6.00	-6.00		
	65	-9.00	+7.00	-6.00		
	70	-7.50	+6.00	-5.00		
	75	-6.00	+5.00	-4.00		
	80	-4.00	+3.50	-3.00		
Transitions® Gen ^S BlueUV Capture™	55	-10.00	+6.00	-6.00		
	60	-10.00	+6.00	-6.00		
	65	-9.00	+6.00	-6.00		
	70	-8.00	+6.00	-5.00		
	75	-7.00	+5.00	-4.00		
EyeProtect System™	55	-10.00	+6.00	-4.00		
	60	-10.00	+6.00	-4.00		
	65	-8.00	+6.00	-4.00		
	70	-7.00	+4.00	-4.00		
	75	-6.00	+5.00	-4.00		
Xperio.	55	-10.00	+6.00	-4.00		
	60	-10.00	+6.00	-4.00		
	65	-8.00	+6.00	-4.00		
	70	-7.00	+4.00	-4.00		

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
XTRACTIVE® NEW GENERATION	Varilux® Comfort Max					
	60	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 5.00 DΔ
	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-8.00	+6.00	-5.00		
	80	-7.00	+5.00	-4.00		
	60	-10.00	+6.00	-6.00		
XTRACTIVE® POLARIZED™	65	-10.00	+6.00	-6.00		
	70	-9.00	+6.00	-6.00		
	75	-8.00	+6.00	-5.00		
	80	-7.00	+5.00	-4.00		

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Varilux® Digitime®



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Prevencia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters	Types	Minimum Depth of Field
Varilux® Digitime™	18mm	Personalised: Screen Distance, Monocular Pupillary Distances, Fitting Heights	Near Mid Room	80cm
Varilux® Digitime™ Fit		Personalised: Screen Distance, Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance		100cm
Varilux® Digitime™ Eyecode		Personalised / Visiooffice only: Screen Distance, Eyecode data (Monocular Pupillary Distances, Fitting Heights & Eye rotation centre distance), Frame Wrap, Pantoscopic Tilt, Back Vertex Distance, Head cape		220cm

Range

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Varilux® Digitime™					
	60	-14.00	+9.00	-6.00	0.75 to 3.50	Combined 6.50 DΔ Standard Fit Eyecode
	65	-12.00	+9.00	-6.00		
	70	-11.00	+9.00	-6.00		
	75	-10.00	+8.00	-5.00		
	80	-9.00	+6.00	-4.00		
1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Varilux® Digitime™					
	65	-12.00	+8.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Standard Fit Eyecode
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
	80	-8.00	+6.00	-4.00		
EyeProtect System™	60	-12.00	+8.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Standard Fit Eyecode
	65	-12.00	+8.00	-6.00		
	70	-10.00	+8.00	-6.00		
	75	-9.00	+8.00	-5.00		
	80	-8.00	+6.00	-4.00		

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Varilux® Digitime™					
	65	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Standard Fit Eyecode
	70	-9.50	+6.00	-6.00		
	75	-8.50	+6.00	-5.00		
	80	-7.25	+5.00	-4.00		
	60	-10.00	+6.00	-6.00		
EyeProtect System™	65	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 6.00 DΔ Standard Fit Eyecode
	70	-9.50	+6.00	-6.00		
	75	-8.50	+6.00	-5.00		
	80	-7.25	+5.00	-4.00		
1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Varilux® Digitime™					
	65	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 5.00 DΔ Standard Fit Eyecode
	70	-9.00	+6.00	-6.00		
	75	-8.00	+6.00	-5.00		
	80	-7.00	+5.00	-4.00		
	60	-10.00	+6.00	-6.00		
EyeProtect System™	65	-10.00	+6.00	-6.00	0.75 to 3.50	Combined 5.00 DΔ Standard Fit Eyecode
	70	-9.00	+6.00	-6.00		
	75	-8.00	+5.00	-5.00		
	80	-7.00	+5.00	-4.00		

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height.

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Essentials

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Prevencia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Essilor Essentials™ Assent	17mm	Monocular Pupillary Distances and Fitting Heights
Essilor Essentials™ Assent Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights

Range

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Essilor Essentials™ Assent					
	65	-12.00	+8.00	+6.00	0.75 to 3.50	Combined 6.50 DA
	70	-12.00	+6.00	+6.00		
	75	-11.00	+4.75	+5.00		
	80	-10.00	+4.75	+4.00		
Transitions Gen [®] S EyeProtect System™	60	-12.00	+9.00	+6.00		
	65	-12.00	+9.00	+6.00		
	70	-11.00	+9.00	+6.00		
	75	-10.00	+8.00	+5.00		
XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED™	60	-14.00	+9.00	+6.00		
	65	-12.00	+9.00	+6.00		
	70	-11.00	+9.00	+6.00		
	75	-10.00	+8.00	+5.00		
Xperio.	65	-12.00	+9.00	+4.00		
	70	-11.00	+9.00	+4.00		
	75	-10.00	+8.00	+4.00		
	80	-9.00	+6.00	+4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Essilor Essentials™ Assent					
	65	-10.00	+6.00	+6.00	0.75 to 3.50	Combined 6.00 DA
	70	-10.00	+4.50	+6.00		
	75	-10.00	+4.50	+5.00		
	80	-9.50	+4.50	+4.00		
EyeProtect System™	60	-12.00	+8.00	+6.00		
	65	-12.00	+8.00	+6.00		
	70	-10.00	+8.00	+6.00		
	75	-9.00	+8.00	+5.00		
Transitions Gen [®] S	60	-12.00	+8.00	+6.00		
	65	-12.00	+8.00	+6.00		
	70	-10.00	+8.00	+6.00		
	75	-9.00	+8.00	+5.00		
XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED™	60	-12.00	+8.00	+6.00		
	65	-12.00	+8.00	+6.00		
	70	-10.00	+8.00	+6.00		
	75	-9.00	+8.00	+5.00		
Xperio.	60	-12.00	+8.00	+4.00		
	65	-10.00	+8.00	+4.00		
	70	-9.00	+8.00	+4.00		
	75	-8.00	+7.00	+4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Essentials

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire™^{HR}

Crizal®
Prevenica®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Essilor Essentials™ Assent	17mm	Monocular Pupillary Distances and Fitting Heights
Essilor Essentials™ Assent Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights

Range

1.59 [*]	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Essilor Essentials™ Assent					
	65	-10.00	+6.00	+6.00	0.75 to 3.50	Combined 6.00 DA
	70	-10.00	+6.00	+5.00		
EyeProtect System™	75	-10.00	+6.00	+4.00		
	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.50	+6.00	+6.00		
Transitions Gen® S	75	-8.50	+6.00	+5.00		
	79	-7.25	+5.00	+4.00		
	65	-10.00	+6.00	+6.00		
	70	-9.50	+6.00	+6.00		
Xperio®	75	-8.50	+6.00	+5.00		
	79	-7.25	+5.00	+4.00		
	60	-10.00	+8.00	+6.00		
	65	-10.00	+7.00	+5.00		
	70	-9.00	+6.00	+4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear	Essilor Essentials™ Assent					
	65	-8.00	+6.00	+6.00	0.75 to 3.50	Combined 5.00 DA
	70	-8.00	+6.00	+6.00		
	75	-8.00	+4.50	+5.00		
EyeProtect System™	80	-5.00	+4.50	+4.00		
	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
Transitions Gen® S	75	-8.00	+6.00	+5.00		
	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
XTRACTIVE® NEW GENERATION	75	-8.00	+6.00	+5.00		
	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
XTRACTIVE® POLARIZED™	75	-8.00	+6.00	+5.00		
	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
Xperio®	75	-8.00	+6.00	+5.00		
	60	-4.00	+6.00	+4.00		
	65	-4.00	+6.00	+4.00		
	70	-4.00	+6.00	+4.00		

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Essentials

Coatings

Crizal[®]
Easy^{Pro}

Crizal[®]
Rock[™]

Crizal[®]
Sapphire^{™HR}

Crizal[®]
Prevenica[®]

Crizal[®]
Sun XProtect[®]

Optifog[®]

Versions

	MFH ⁽³⁾	Parameters
Essilor Essentials™ Accord	17mm	Monocular Pupillary Distances and Fitting Heights
Essilor Essentials™ Accord Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights

Range

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Essilor Essentials™ Accord						
Clear Transitions Gen[®]S EyeProtect System[™]	60	-12.00	+9.00	+6.00	0.75 to 3.50	Combined 6.50 DΔ Regular Short
	65	-12.00	+9.00	+6.00		
	70	-11.00	+9.00	+6.00		
	75	-10.00	+8.00	+5.00		
	80	-9.00	+6.00	+4.00		
XTRACTIVE[®]* NEW GENERATION XTRACTIVE[®] POLARIZED*	60	-14.00	+9.00	+6.00		
	65	-12.00	+9.00	+6.00		
	70	-11.00	+9.00	+6.00		
	75	-10.00	+8.00	+5.00		
	80	-9.00	+6.00	+4.00		
Xperio. • • •	65	-12.00	+9.00	+4.00		
	70	-11.00	+9.00	+4.00		
	75	-10.00	+8.00	+4.00		
	80	-9.00	+6.00	+4.00		

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Essilor Essentials™ Accord						
Clear Transitions Gen[®]S EyeProtect System[™]	60	-12.00	+8.00	+6.00	0.75 to 3.50	Combined 6.00 DΔ Regular Short
	65	-12.00	+8.00	+6.00		
	70	-10.00	+8.00	+6.00		
	75	-9.00	+8.00	+5.00		
	80	-8.00	+6.00	+4.00		
XTRACTIVE[®]* NEW GENERATION XTRACTIVE[®] POLARIZED*	60	-12.00	+8.00	+6.00		
	65	-12.00	+8.00	+6.00		
	70	-10.00	+8.00	+6.00		
	75	-9.00	+8.00	+5.00		
	80	-8.00	+6.00	+4.00		
Xperio. • • •	65	-12.00	+8.00	+4.00		
	70	-10.00	+8.00	+4.00		
	75	-9.00	+8.00	+4.00		
	80	-8.00	+7.00	+4.00		

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Essentials

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Prevenia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Essilor Essentials™ Accord	17mm	Monocular Pupillary Distances and Fitting Heights
Essilor Essentials™ Accord Short	14mm	Short corridor: Addition up to 3.00 D. Monocular Pupillary Distances and Fitting Heights

Range

1.59*	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear EyeProtect System..	60	-10.00	+6.00	+6.00	0.75 to 3.50	Combined 6.50 DΔ
	65	-10.00	+6.00	+6.00		
	70	-9.50	+6.00	+6.00		
	75	-8.50	+6.00	+5.00		
	80	-7.25	+5.00	+4.00		
Transitions Gen [®] S™	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.50	+6.00	+6.00		
	75	-8.50	+6.00	+5.00		
	79	-7.25	+5.00	+4.00		
Xperio®	65	-10.00	+6.00	+6.00		
	70	-9.50	+6.00	+6.00		
	75	-8.50	+6.00	+5.00		
	80	-7.25	+5.00	+4.00		

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽⁴⁾	Versions
Clear Transitions Gen [®] S™	60	-10.00	+6.00	+6.00	0.75 to 3.50	Combined 5.00 DΔ
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
	75	-8.00	+6.00	+5.00		
	80	-7.00	+5.00	+4.00		
XTRACTIVE® NEW GENERATION	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
	75	-8.00	+6.00	+5.00		
	80	-7.00	+5.00	+4.00		
XTRACTIVE® POLARIZED™	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
	75	-8.00	+6.00	+5.00		
EyeProtect System..	60	-10.00	+6.00	+6.00		
	65	-10.00	+6.00	+6.00		
	70	-9.00	+6.00	+6.00		
	75	-8.00	+6.00	+5.00		
Xperio®	65	-10.00	+6.00	+4.00		
	70	-8.00	+6.00	+4.00		
	75	-7.00	+4.00	+4.00		

*Note: Transitions XTRACTIVE New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Eyezen™ Boost



Coatings

**Crizal®
EasyPro**

**Crizal®
Rock™**

**Crizal®
Sapphire™^{HR}**

**Crizal®
Previncia®**

**Crizal®⁽²⁾
Sun XProtect**

Optifog®

Versions

Eyezen™ Boost Rx is available in:	Regular: Monocular Pupillary Distances and Fitting Heights
	Fit: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
	Eyecode: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

	MFH ⁽³⁾	Shift
Eyezen™ Boost Plus	17mm	0.40 D
Eyezen™ Boost Pro	17mm	0.60 D
Eyezen™ Boost Max	17mm	0.85 D
Eyezen™ Boost Expert	17mm	1.10 D

Range

1.74	Ø	Range	Cyl Max	Prism	Versions
Clear Transitions Gen®S EyeProtect System...	Eyezen™ Boost				
	60	-18.00	+13.00	-6.00	Combined 6.50 ΔΔ Plus Pro Max
	65	-15.00	+13.00	-6.00	
	70	-12.00	+9.00	-6.00	
	75	-10.00	+8.00	-5.00	
	80	-9.00	+6.00	-4.00	
1.67	Ø	Range	Cyl Max	Prism	Versions
Clear Transitions Gen®S EyeProtect System...	Eyezen™ Boost				
	65	-12.00	+9.00	-6.00	Combined 6.50 ΔΔ Plus Pro Max
	70	-11.00	+9.00	-6.00	
	75	-10.00	+8.00	-5.00	
	80	-9.00	+6.00	-4.00	
	80	-9.00	+6.00	-4.00	
XTRACTIVE® NEW GENERATION EyeProtect System...	60	-14.00	+9.00	-6.00	Combined 6.50 ΔΔ Plus Pro Max
	65	-12.00	+9.00	-6.00	
	70	-11.00	+9.00	-6.00	
	75	-10.00	+8.00	-5.00	
	80	-9.00	+6.00	-4.00	

1.60	Ø	Range	Cyl Max	Prism	Versions
Clear Transitions Gen®S EyeProtect System...	Eyezen™ Boost				
	65	-12.00	+8.00	-6.00	Combined 6.00 ΔΔ Plus Pro Max
	70	-10.00	+8.00	-6.00	
	75	-9.00	+8.00	-5.00	
	80	-8.00	+6.00	-4.00	
	80	-8.00	+6.00	-4.00	
XTRACTIVE® NEW GENERATION EyeProtect System...	60	-12.00	+8.00	-6.00	Combined 6.00 ΔΔ Plus Pro Max
	65	-12.00	+8.00	-6.00	
	70	-10.00	+8.00	-6.00	
	75	-9.00	+8.00	-5.00	
	80	-8.00	+6.00	-4.00	
1.59	Ø	Range	Cyl Max	Prism	Versions
Clear Transitions Gen®S EyeProtect System...	Eyezen™ Boost				
	60	-10.00	+6.00	-6.00	Combined 6.00 ΔΔ Plus Pro Max
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-8.50	+6.00	-5.00	
	75	-8.50	+6.00	-5.00	

*Note: Transitions XTRACTIVE New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.50	Ø	Range	Cyl Max	Prism	Versions
Clear Transitions Gen®S EyeProtect System...	Eyezen™ Boost				
	60	-10.00	+6.00	-6.00	Combined 5.00 ΔΔ Plus Pro Max
	65	-10.00	+6.00	-6.00	
	70	-9.00	+6.00	-6.00	
	75	-8.00	+6.00	-5.00	
	80	-7.00	+5.00	-4.00	
XTRACTIVE® NEW GENERATION EyeProtect System...	60	-10.00	+6.00	-6.00	Combined 5.00 ΔΔ Plus Pro Max
	65	-10.00	+6.00	-6.00	
	70	-9.00	+6.00	-6.00	
	75	-8.00	+6.00	-5.00	
	75	-8.00	+6.00	-5.00	

(1) Crizal® Easy Pro coating is not available in 1.74 and 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height



Eyezen™ Start

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire™^{HR}

Crizal®
Prevencia®

Crizal®⁽²⁾
Sun XProtect

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Eyezen™ Start	N / A	Monocular Pupillary Distances and Fitting Heights
Eyezen™ Start Fit		Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Eyezen™ Start Eyezen		Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.74	Ø	Range	Cyl Max	Prism ⁽⁴⁾	Versions
Clear	Eyezen™ Start				
	60	-18.00	13.00	6.00	Combined 6.50 DΔ
	65	-18.00	13.00	6.00	
	70	-11.00	9.00	6.00	
	75	-11.00	9.00	6.00	
	80	-11.00	6.00	6.00	
EyeProtect System™	60	-18.00	13.00	6.00	
	65	-18.00	13.00	6.00	
	70	-11.00	9.00	6.00	
	75	-11.00	9.00	6.00	
	80	-11.00	9.00	6.00	
Transitions® Gen®S	55	-12.00	6.00	6.00	
	60	-12.00	6.00	6.00	
	65	-10.00	6.00	6.00	
	70	-9.00	6.00	6.00	
	75	-8.00	6.00	6.00	
	80	-7.00	4.00	6.00	

1.67	Ø	Range	Cyl Max	Prism ⁽⁴⁾	Versions
Clear BlueUV Capture™	Eyezen™ Start				
	60	-14.00	+9.00	-6.00	Combined 6.50 DΔ
	65	-14.00	+9.00	-6.00	
	70	-13.00	+9.00	-6.00	
	75	-12.00	+9.00	-6.00	
	80	-11.00	+8.00	-6.00	
EyeProtect System™ Xperio®	65	-12.00	+9.00	-4.00	
	70	-11.00	+9.00	-4.00	
	75	-10.00	+8.00	-4.00	
	80	-9.00	+6.00	-4.00	
Transitions® gen®S XTRACTIVE®** NEW GENERATION XTRACTIVE® POLARIZED	60	-14.00	+9.00	-6.00	
	65	-14.00	+9.00	-6.00	
	70	-13.00	+9.00	-6.00	
	75	-12.00	+9.00	-6.00	
	80	-11.00	+8.00	-6.00	

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal® Easy Pro coating is not available in 1.74 and 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

(4) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Eyezen™ Start

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®⁽²⁾
Sun XProtect

Optifog®

Versions


	MFH ⁽³⁾	Parameters
Eyezen™ Start	N / A	Monocular Pupillary Distances and Fitting Heights
Eyezen™ Start Fit		Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Eyezen™ Start Eyecode		Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.60	Ø	Range	Cyl Max	Prism ⁽⁴⁾	Versions
Clear BlueUV Capture™	Eyezen™ Start				
	60	-12.00	+8.00	-6.00	Combined 6.00 DΔ
	65	-12.00	+8.00	-6.00	
	70	-11.00	+8.00	-6.00	
	75	-10.50	+8.00	-6.00	
	80	-10.00	+6.00	-5.00	
EyeProtect System™ Xperio®	60	-12.00	+8.00	-6.00	
	65	-12.00	+8.00	-6.00	
	70	-11.00	+8.00	-6.00	
	75	-10.50	+8.00	-6.00	
Transitions Gen®S XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED™	60	-12.00	+8.00	-6.00	
	65	-12.00	+8.00	-6.00	
	70	-11.00	+8.00	-6.00	
	75	-10.50	+8.00	-6.00	
	80	-10.00	+6.00	-5.00	
					Regular Fit Eyecode

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.59	Ø	Range	Cyl Max	Prism ⁽⁴⁾	Versions
Clear	Eyezen™ Start				
	60	-10.00	+6.00	-6.00	Combined 6.00 DΔ
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-9.00	+6.00	-6.00	
	80	-8.50	+5.50	-6.00	
BlueUV Capture™	60	-10.00	+6.00	-6.00	
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-9.00	+6.00	-6.00	
EyeProtect System™ Xperio®	60	-10.00	+6.00	-6.00	
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-9.00	+6.00	-6.00	
Transitions Gen®S	60	-10.00	+6.00	-6.00	
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-9.00	+6.00	-6.00	
	80	-8.50	+5.50	-6.00	
					Regular Fit Eyecode



Eyezen™ Start

Coatings

**Crizal®
EasyPro**⁽¹⁾

**Crizal®
Rock™**

**Crizal®
Sapphire™^{HR}**

**Crizal®
Previncia®**

**Crizal®
Sun XProtect**⁽²⁾

Optifog®

Versions

	MFH ⁽³⁾	Parameters
Eyezen™ Start	N / A	Monocular Pupillary Distances and Fitting Heights
Eyezen™ Start Fit		Personalised: Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance
Eyezen™ Start Eyecode		Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance & Eyecode (Eye Rotation Centre distance and Head Cape)

Range

1.50	Ø	Range	Cyl Max	Prism ⁽⁴⁾	Versions
Clear BlueUV Capture™	Eyezen™ Start				
	60	-10.00	+6.00	-6.00	Combined 5.00 DΔ
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-8.50	+6.00	-5.00	
	80	-7.50	+6.00	-5.00	
EyeProtect System™ Xperio®	60	-10.00	+6.00	-6.00	
	65	-10.00	+6.00	-6.00	
	70	-9.00	+6.00	-6.00	
	75	-8.50	+6.00	-5.00	
Transitions Gen#S™ XTRACTIVE® NEW GENERATION XTRACTIVE® POLARIZED	60	-10.00	+6.00	-6.00	
	65	-10.00	+6.00	-6.00	
	70	-9.50	+6.00	-6.00	
	75	-8.50	+6.00	-5.00	

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.



Coatings

Crizal[®]
Easy^{Pro}
Crizal[®]
Rock[™]
Crizal[®]
Sapphire^{HR}
Crizal[®]
Previncia[®]
Crizal[®]
Sun XProtect⁽²⁾
Optifog[®]

Versions

	MFH ⁽³⁾	Parameters
SV 360 Fit Aspheric	N / A	Monocular Pupillary Distances, Fitting Heights, Frame Wrap, Pantoscopic Tilt & Back Vertex Distance

Range

1.74	Ø	Range		Cyl Max
Clear	Essilor SV 360 Fit			
	60	-14.00	+10.00	-6.00
	65	-14.00	+10.00	-6.00
	70	-12.00	+6.00	-6.00
	75	-11.00	+5.00	-6.00
EyeProtect System	55	-14.00	+10.00	-6.00
	60	-14.00	+10.00	-6.00
	65	-16.00	+10.00	-6.00
	70	-12.00	+10.00	-6.00
	75	-11.50	+8.00	-5.00
	80	-11.00	+7.00	-4.00

1.67	Ø	Range	Cyl Max	
Clear	Essilor SV 360 Fit			
	65	-13.00	+8.00	-6.00
	70	-11.00	+7.00	-6.00
	75	-8.00	+4.00	-6.00
EyeProtect System™	55/60/65	-13.00	+9.00	-6.00
	70	-11.00	+9.00	-6.00
	75	-10.00	+8.50	-5.00
	80	-9.00	+7.50	-4.00
	55/60/65	-13.00	+9.00	-6.00
	70	-11.00	+9.00	-6.00
	75	-10.00	+8.50	-5.00
	80	-9.00	+7.50	-4.00

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal[®] Easy Pro coating is not available in 1.74 and 1.67 index.

(2) Crizal[®] Sun XProtect[™] is only available on tinted and polarised lenses (not with mirror coatings).

(3) Minimum Fitting Height

1.60	Ø	Range	Cyl Max	
Clear	Essilor SV 360 Fit			
	65	-13.00	+8.00	-6.00
	70	-11.00	+5.00	-6.00
	75	-8.00	+4.00	-6.00
EyeProtect System...	55/60	-13.00	+8.00	-6.00
	65	-11.50	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-9.00	+8.00	-5.00
	80	-8.00	+6.00	-4.00
Transitions gen[®]S 	55/60	-13.00	+8.00	-6.00
	65	-11.50	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-9.00	+8.00	-5.00
	80	-8.00	+6.00	-4.00

1.59	Ø	Range	Cyl Max	
Clear EyeProtect System..	Essilor SV 360 Fit			
	60	-10.00	+8.00	-6.00
	65	-10.00	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-10.00	+8.00	-6.00
Transitions Gen[®]S 	60	-10.00	+8.00	-6.00
	65	-10.00	+8.00	-6.00
	70	-10.00	+8.00	-6.00

1.50	Ø	Range	Cyl Max	
Clear EyeProtect System..	Essilor SV 360 Fit			
	60	-12.00	+10.00	-6.00
	65	-10.00	+10.00	-5.00
	70	-7.00	+8.00	-4.00
	75	-3.50	+6.00	-4.00
Transitions Gen[®]S 	60	-11.00	+10.00	-6.00
	70	-9.00	+6.00	-6.00
	55	-12.00	+10.00	-6.00
	65	-10.00	+10.00	-6.00
	75	-8.00	+4.50	-5.00



Coatings

Crizal[®]
Easy^{Pro}

Crizal[®]
Rock[™]

Crizal[®]
Sapphire^{HR}

Crizal[®]
Previncia[®]

Crizal[®]
Sun XProtect⁽²⁾

Optifog[®]

Range					
1.74	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Aspheric				Combined 6.50 DΔ
	55/60	-18.00	+10.00	-6.00	
	65	-16.00	+10.00	-6.00	
	70	-13.00	+6.00	-6.00	
	75	-10.00	+2.00	-6.00	
	80	-8.00	+1.00	-6.00	
Transitions Gen [®] S [™]	55	-18.00	+10.00	-6.00	
	60	-18.00	+10.00	-6.00	
	65	-16.00	+10.00	-6.00	
	70	-12.00	+6.00	-6.00	
	75	-10.00	+2.00	-6.00	
	80	-8.00	+2.00	-6.00	
1.60	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Aspheric				Combined 6.00 DΔ
	55/65	-15.00	+10.50	+10.50	
	70	-15.00	+10.50	+10.50	
	75	-15.00	+10.50	+10.50	
Transitions Gen [®] S [™]	55	-13.00	+8.00	-6.00	
	60	-13.00	+8.00	-6.00	
	65	-11.50	+8.00	-6.00	
	70	-10.00	+8.00	-6.00	
	75	-9.00	+8.00	-5.00	
	80	-8.00	+6.00	-4.00	

Out of Range:

Clear and Transitions available up to +10.50 to combined -15.00
/ Cyl -10.50 Combined 10.50 DΔ

1.67	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Aspheric				Combined 6.50 DΔ
	65	-17.00	+8.50	-13.00	
	70	-17.00	+8.50	-13.00	
	75	-17.00	+8.50	-13.00	
	80	-17.00	+8.50	-13.00	
Transitions Gen [®] S [™]	55	-13.00	-8.00	-6.00	
	60	-13.00	-8.50	-6.00	
	65	-13.00	-8.50	-6.00	
	70	-11.00	-4.00	-6.00	
	75	-10.00	-3.00	-6.00	
	80	-9.00	-3.00	-4.00	

1.59	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Aspheric				Combined 6.00 DΔ
	75	-15.00	+9.00	-6.00	
Transitions Gen [®] S [™]	55	-10.00	+8.00	-6.00	
	60	-10.00	+8.00	-6.00	
	65	-10.00	+8.00	-6.00	
	70	-9.00	+8.00	-5.00	
	74	-8.50	+6.00	-5.00	
	79	-8.00	+4.50	-4.00	

(1) Crizal[®] Easy Pro coating is not available in 1.67 index.

(2) Crizal[®] Sun XProtect[™] is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Coatings

Crizal⁽¹⁾
Easy^{Pro}
Crizal[®]
Rock[™]
Crizal[®]
Sapphire^{™HR}
Crizal[®]
Previncia[®]
Crizal^{®(2)}
Sun XProtect

Optifog[®]

Range

1.67	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Spheric				
	65	-15.00	+10.00	-10.00	Combined 6.50 DΔ
	70	-15.00	+10.00	-10.00	
	75	-15.00	+10.00	-10.00	
	80	-15.00	+10.00	-10.00	
EyeProtect System	55	-13.00	+9.00	-6.00	
	60	-13.00	+9.00	-6.00	
	65	-13.00	+9.00	-6.00	
	70	-11.00	+9.00	-6.00	
	75	-10.00	+8.50	-5.00	
	80	-9.00	+7.50	-4.00	
Transitions Gen [®] S [™]	55/60	-13.00	+9.00	-6.00	
	65	-13.00	+9.00	-6.00	
	70	-11.00	+9.00	-6.00	
	75	-10.00	+8.50	-5.00	
	80	-9.00	+7.50	-4.00	
XTRACTIVE [®] * NEW GENERATION  XTRACTIVE [®] * POLARIZED [®] 	55	-13.00	+9.00	-6.00	
	60	-13.00	+9.00	-6.00	
	65	-13.00	+9.00	-6.00	
	70	-11.00	+9.00	-6.00	
	75	-10.00	+8.50	-5.00	
	80	-9.00	+7.50	-4.00	
 Xperio	60	-13.00	+8.00	-4.00	
	65	-13.00	+8.00	-4.00	
	70	-11.00	+8.00	-4.00	
	75	-10.00	+8.00	-4.00	
	80	-7.50	+7.00	-4.00	

1.60	Ø	Range	Cyl Max	Prism ⁽³⁾
Clear	Spheric			
	65	-13.00	+10.00	-10.00
	70	-13.00	+10.00	-10.00
	75	-13.00	+10.00	-10.00
	80	-13.00	+10.00	-10.00
EyeProtect System..	55	-13.00	+8.00	-6.00
	60	-13.00	+8.00	-6.00
	65	-11.50	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-9.00	+8.00	-5.00
	80	-8.00	+6.00	-4.00
 Transitions Gen [®] S	55	-13.00	+8.00	-6.00
	60	-13.00	+8.00	-6.00
	65	-11.50	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-9.00	+8.00	-5.00
	80	-8.00	+6.00	-4.00
 XTRACTIVE [®] * NEW GENERATION  XTRACTIVE [®] * POLARIZED [®]	55	-13.00	+8.00	-6.00
	60	-13.00	+8.00	-6.00
	65	-11.50	+8.00	-6.00
	70	-10.00	+8.00	-6.00
	75	-9.00	+8.00	-5.00
	80	-8.00	+6.00	-4.00
	65	-12.00	+8.00	-4.00
	70	-10.00	+8.00	-4.00
	75	-9.00	+8.00	-4.00
	80	-8.00	+7.00	-4.00

***Note:** Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

Out of Range: Transitions Grey available up to +15.00 to combined -15.00 / Cyl -10.50 combined 10.00 DΔ
Xperio[®] available up to +8.00 to combined -15.00 / Cyl -10.50 combined 10.00 DΔ

(1) Crizal[®] Easy Pro coating is not available in 1.67 index.

(2) Crizal[®] Sun XProtect[™] is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details. Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter. Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.



Coatings

Crizal[®]
Easy^{Pro}
Crizal[®]
Rock[™]
Crizal[®]
Sapphire^{HR}
Crizal[®]
Previncia[®]
Crizal[®]
Sun XProtect⁽²⁾
Optifog[®]

Range					
1.59*	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Spheric				Combined 6.00 Δ
	55/60	-15.00	+14.50	-10.50	
	65/70	-15.00	+14.50	-10.50	
	75	-15.00	+14.50	-10.50	
EyeProtect System..	55	-10.00	+8.00	-6.00	
	60	-10.00	+8.00	-6.00	
	65	-10.00	+8.00	-6.00	
	70	-9.00	+8.00	-5.00	
	74	-8.50	+6.00	-5.00	
	79	-8.00	+4.50	-4.00	
Transitions Gen#S	55	-10.00	+8.00	-6.00	
	60	-10.00	+8.00	-6.00	
	65	-10.00	+8.00	-6.00	
	70	-9.00	+8.00	-6.00	
	74	-8.50	+6.00	-5.00	
	79	-8.00	+4.50	-4.00	
Xperio.	55/60	-15.00	+9.50	+10.50	
	65/70	-15.00	+9.50	+10.50	
	75/76	-15.00	+9.50	+10.50	
	60	-10.00	+6.00	-4.00	
XP TINT	65	-10.00	+6.00	-4.00	
	70	-10.00	+6.00	-4.00	
	74	-4.00	+4.00	-4.00	

1.50	Ø	Range	Cyl Max	Prism ⁽³⁾	
Clear	Spheric				Combined 5.00 Δ
	70	-12.00	+9.75	-10.00	
	76	-12.00	+9.75	-10.00	
EyeProtect System..	55	-12.00	+10.00	-6.00	
	60	-11.00	+10.00	-6.00	
	65	-10.00	+10.00	-6.00	
	70	-9.00	+6.00	-6.00	
	75	-8.00	+4.50	-6.00	
Transitions Gen#S	55	-12.00	+10.00	-6.00	
	60	-11.00	+10.00	-6.00	
	65	-10.00	+10.00	-6.00	
	70	-9.00	+6.00	-6.00	
XTRACTIVE® NEW GENERATION	75	-8.00	+4.50	-6.00	
	55	-12.00	+10.00	-6.00	
	60	-11.00	+10.00	-6.00	
	65	-10.00	+10.00	-6.00	
XTRACTIVE® POLARIZED®	70	-9.00	+6.00	-6.00	
	75	-8.00	+5.00	-5.00	
	Xperio.	80	-8.00	+7.00	

***Note:** Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

(1) Crizal[®] Easy Pro coating is not available in 1.67 index.

(2) Crizal[®] Sun XProtect[™] is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details. Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter. Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.



S O L A R F L A I R

Coatings

S O L A R F L A I R
Sun XProtect

Versions

Base Curve

5, 6, 7 and 8. Minimum Fitting Height: 17mm

Range

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair Progressive					
	80	-6.00	+6.00	-4.00	0.75 to 3.50	5 / 6
	85	-5.00	+5.00	-3.00	Combined 6.50 DΔ	5 / 6 / 7 / 8
	90	-4.00	+4.00	-3.00		
	95	-3.00	+3.00	-2.00		

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair Progressive					
	80	-4.00	+4.00	-4.00	0.75 to 3.50	5 / 6
	90	-4.00	+4.00	-3.00	Combined 6.00 DΔ	5 / 6 / 7 / 8
	95	-3.00	+3.00	-2.00		

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair Progressive					
	80	-5.00	+5.00	-4.00	0.75 to 3.50	5 / 6
	85	-4.00	+4.00	-3.00	Combined 6.00 DΔ	5 / 6 / 7 / 8
	90	-3.00	+3.00	-3.00		
	95	-2.00	+3.00	-2.00		

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair Progressive					
	80	-4.00	+3.00	-4.00	Combined 5.00 DΔ	5 / 6
	85	-3.50	+3.00	-3.00		5 / 6 / 7 / 8
	90	-3.00	+3.00	-3.00		
	95	-2.00	+2.00	-2.00		7 / 8

Tints



Grey



Brown



Green

Grey Blue⁽²⁾

Mirrors



Gold Pink



Silver



Bronze



Silver Shadow



Violet



Turquoise Blue



Forest Green



Ocean Blue

(1) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

(2) Not available in polarised.



S O L A R F L A I R

Coatings

S O L A R F L A I R
Sun XProtect

Versions

Base Curve

5, 6, 7 and 8.

Range

1.67	Ø	Range	Cyl Max	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair SV				
	80	-6.00	+6.00	-4.00	5 / 6
	85	-5.00	+5.00	-3.00	
	90	-4.00	+4.00	-3.00	5 / 6 / 7 / 8
	95	-3.00	+3.00	-2.00	7 / 8
Combined 6.50 DΔ					

1.59	Ø	Range	Cyl Max	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair SV				
	80	-4.00	+4.00	-4.00	5 / 6
	90	-4.00	+4.00	-3.00	
	95	-3.00	+3.00	-2.00	7
Combined 6.00 DΔ					

1.60	Ø	Range	Cyl Max	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair SV				
	80	-5.00	+5.00	-4.00	5 / 6
	85	-4.00	+4.00	-3.00	
	90	-3.00	+3.00	-3.00	5 / 6 / 7 / 8
	95	-2.00	+3.00	-2.00	7 / 8
Combined 6.00 DΔ					

1.50	Ø	Range	Cyl Max	Prism ⁽¹⁾	Base Curve
TINT (Polar / Non-Polar) MIRROR (Polar / Non-Polar)	Solarflair SV				
	80	-4.00	+3.00	-4.00	5 / 6
	85	-3.50	+3.00	-3.00	
	90	-3.00	+3.00	-3.00	5 / 6 / 7 / 8
	95	-2.00	+2.00	-2.00	7
Combined 5.00 DΔ					

Tints



Grey



Brown



Green

Grey Blue⁽²⁾

Mirrors



Gold Pink



Silver



Bronze



Silver Shadow



Violet



Turquoise Blue



Forest Green



Ocean Blue

(1) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

(2) Not available in polarised.

Stellest™



Coatings

Crizal®
Rock™

Versions

	MFH ⁽¹⁾	Parameters
Stellest	N / A	Monocular Pupillary Distances and Fitting Heights


Range

1.59	Ø	Range	Prism ⁽²⁾	Cyl Max	
Clear	Stellest				
	55	-12.00	+2.00	Combined 5.00 Δ	-4.00
	60	-11.00	+2.00		
	65	-10.00	+2.00		
	70	-9.00	+2.00		

Essilor® Stellest® in Clear lenses are available with Crizal® Rock coatings.

NEW - Essilor® Stellest® with Sun Tints

Essilor® Stellest® with Sun Tints is available in 65mm and 70mm diameter;
Sph -6.00 to 0.00, Cyl -2.00 with Crizal® Sun XProtect.



Note: Category 3 colours are not compliant with ISO standards for road use and driving at night. Emerald Blue tint colour is not compliant for driving during daylight as well.

(1) Minimum Fitting Height

(2) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

myopilux®



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Range

Myopilux® Pro

1.67	Ø	Range	Cyl Max	ADD Min	ADD Max
Clear	Myopilux® Pro				
	60	-8.00	+0.00	+4.00	1.50
	65	-8.00	+0.00		
	70	-8.00	+0.00		
	75	-8.00	+0.00		
	80	-8.00	+0.00		

1.60	Ø	Range	Cyl Max	ADD Min	ADD Max
Clear	Myopilux® Pro				
	60	-8.00	+0.00	+4.00	1.50
	65	-8.00	+0.00		
	70	-8.00	+0.00		
	75	-8.00	+0.00		
	80	-8.00	+0.00		

1.59	Ø	Range	Cyl Max	ADD Min	ADD Max
Clear	Myopilux® Pro				
	60	-8.00	+0.00	+4.00	1.50
	65	-8.00	+0.00		
	70	-8.00	+0.00		
	75	-8.00	+0.00		
	80	-7.25	+0.00		

1.56	Ø	Range	Cyl Max	ADD Min	ADD Max
Clear	Myopilux® Pro				
	55	-8.00	+0.00	+4.00	1.50
	60	-8.00	+0.00	+4.00	
	65	-8.00	+0.00	+3.00	
	70	-7.00	+0.00	+2.00	

1.50	Ø	Range	Cyl Max	ADD Min	ADD Max
Clear	Myopilux® Pro				
	60	-8.00	+0.00	+4.00	1.50
	65	-8.00	+0.00	+4.00	
	70	-8.00	+0.00	+4.00	
	75	-8.00	+0.00	+2.00	
	80	-7.00	+0.00	+1.00	

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

Varilux® Stylistic™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

Base Curve	5, 6, 7 and 8. Minimum Fitting Height: 17mm
Parameters	Measurements Required: Wrap measurement required. Tilt is optional (highly recommended)

Range

1.67	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾	Base Curve	
Clear Transitions Gen[®]S™ Xperio.	Varilux® Stylistic™						
	80	-6.00	+6.00	-4.00	0.75 to 3.50	Combined 6.50 DΔ	5 / 6
	85	-5.00	+5.00	-3.00	0.75 to 3.00		5 / 6 / 7 / 8
	90	-4.00	+4.00	-3.00			7 / 8
	95	-3.00	+3.00	-2.00			
XTRACTIVE[®] NEW GENERATION	80	-6.00	+6.00	-4.00	0.75 to 3.50	5 / 6	

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.60	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾	Base Curve	
Clear Xperio. 	Varilux® Stylistic™						
	80	-5.00	+5.00	-4.00	0.75 to 3.50	Combined 6.00 DΔ	5 / 6
	85	-4.00	+4.00	-3.00	0.75 to 3.00		
	90	-3.00	+3.00	-3.00			5 / 6 / 7 / 8
	95	-2.00	+3.00	-2.00			7 / 8
XTRACTIVE® NEW GENERATION 	80	-5.00	+5.00	+4.00	0.75 to 3.50	5 / 6	

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾	Base Curve
Clear Transitions Gen#S™ Xperio.	80	-4.00 +4.00	-4.00	0.75 to 3.50	Combined 6.00 DΔ	5 / 6
	85	-4.00 +4.00	-4.00			
	90	-4.00 +4.00	-3.00			5 / 6
	95	-4.00 +4.00	-3.00			5 / 6 / 7 / 8

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾	Base Curve
Clear Transitions Gen#S™ Xperio.	80	-4.00 +4.00	-4.00	0.75 to 3.50	Combined 5.00 DΔ	5 / 6
	85	-4.00 +4.00	-3.00			
	90	-4.00 +4.00	-3.00			5 / 6 / 7 / 8
	95	-4.00 +4.00	-3.00			5 / 6
XTRACTIVE® NEW GENERATION	76	-4.00 +4.00	-4.00	0.75 to 3.50	Combined 5.00 DΔ	5 / 6
	80	-4.00 +3.00	-4.00			
	85	-3.50 +3.00	-3.00			5 / 6
	90	-3.00 +3.00	-3.00			5 / 6 / 7 / 8
Xperio.	95	+2.00 +2.00	-2.00			7 / 8

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

Stylistic™



Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

Base Curve	5, 6, 7 and 8.
Parameters	Measurements Required: Wrap measurement required. Tilt is optional (highly recommended)

Range

1.67	Ø	Range	Cyl Max	Prism ⁽³⁾	Base Curve	
Clear Transitions Gen [®] S [™]  Xperio. 	Single Vision Stylistic™					
	80	-6.00	+6.00	-4.00	Combined 6.50 Δ	5 / 6
	85	-5.00	+5.00	-3.00		
	90	-4.00	+4.00	-3.00		
	95	-3.00	+3.00	-2.00		
	XTRACTIVE [®] NEW GENERATION 	80	-6.00	+6.00	-4.00	5 / 6

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.60	Ø	Range	Cyl Max	Prism ⁽⁵⁾	Base Curve	
Clear Xperio® 	Single Vision Stylistic™					
	80	-5.00	+5.00	-4.00	Combined 6.00 DΔ	5 / 6
	85	-4.00	+4.00	-3.00		
	90	-3.00	+3.00	-3.00		5 / 6 / 7 / 8
	95	-2.00	+3.00	-2.00		7 / 8
XTRACTIVE®* NEW GENERATION 	80	-5.00	+5.00	-4.00		5 / 6

*Note: Transitions XTRActive New Generation Lenses with Mirror coatings are not suitable for night driving. See **Transitions page** for available colours.

1.59	Ø	Range	Cyl Max	Prism ⁽³⁾	Base Curve	
Clear Transitions Gen[®]S™ 	Single Vision Stylistic™					
	80	-4.00	+4.00	-4.00	Combined 6.00 DΔ	5 / 6
	85	-4.00	+4.00	-3.00		
	90	-4.00	+4.00	-3.00		
	95	-2.00	+2.00	-2.00		7
Xperio. 	80	-4.00	+4.00	-4.00	Combined 6.00 DΔ	
	85	-4.00	+4.00	-3.00		
	90	-4.00	+4.00	-3.00		
	95	-3.00	+3.00	-2.00		7

1.50	Ø	Range	Cyl Max	Prism ⁽³⁾	Base Curve
Clear	Single Vision Stylistic™				
	80	-4.00	+4.00	-4.00	5 / 6
	80	-4.00	+4.00	-4.00	5 / 6
	85	-4.00	+4.00	-3.00	
	90	-4.00	+4.00	-3.00	
XTRACTIVE® NEW GENERATION 	76	-4.00	+4.00	-4.00	5 / 6
	80	-4.00	+3.00	-4.00	5 / 6
	85	-3.50	+3.00	-3.00	5 / 6 / 7 / 8
	90	-3.00	+3.00	-3.00	
	95	+2.00	+2.00	-2.00	

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Interview®

Coatings

Crizal®
Easy^{Pro}

Crizal®
Rock™

Crizal®
Sapphire^{HR}

Crizal®
Previncia®

Crizal®
Sun XProtect⁽²⁾

Optifog®

Versions

Interview Low	0.80D Shift
Interview High	1.30D Shift

Range

1.67	Ø	Range		Prism ⁽³⁾	Cyl Max	Versions
Clear	Interview					
	75	-7.00	+7.00	Combined 6.50 ΔΔ	-6.00	Low, High

Out of range: 8.00 -9.00 / cyl -12.00

1.60	Ø	Range		Prism ⁽³⁾	Cyl Max	Versions
Clear	Interview					
	80	-6.00	+5.00	Combined 6.00 ΔΔ	-4.00	Low, High

Out of range: 7.00 -7.00 / cyl -10.50

1.59	Ø	Range		Prism ⁽³⁾	Cyl Max	Versions
Clear	Interview					
	80	-6.00	+5.00	Combined 6.00 DΔ	-4.00	Low, High

Out of range: 6.50 -7.00 / cyl -10.50

1.50	Ø	Range		Prism ⁽³⁾	Cyl Max	Versions
Clear	Interview					
	80	-6.00	+5.00	Combined 5.00 ΔΔ	-4.00	Low, High

(1) Crizal® Easy Pro coating is not available in 1.67 index.

(2) Crizal® Sun XProtect™ is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.



Coatings

Crizal[®]
Easy^{Pro}

Crizal[®]
Rock[™]

Crizal[®]
Sapphire^{™HR}

Crizal[®]
Previncia[®]

Crizal[®]
Sun XProtect^{®(2)}

Optifog[®]

BiFocal Range

Flat Top 28

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Bifocal				
	60	-12.00	+9.00	-6.00	1.00 to 3.00
	65	-12.00	+8.00	-6.00	
	70	-12.00	+7.00	-6.00	
	75	-12.00	+5.50	-6.00	
Combined 4.50 DΔ					
Transitions [®] Signature GEN8	75	-10.00	+7.00	-4.00	1.00 to 3.50

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Bifocal				
	60	-10.50	+9.50	-6.00	0.75 to 4.00
	65	-10.50	+9.50	-6.00	
	70	-10.50	+9.50	-6.00	
	75	-10.50	+9.50	-6.00	
Combined 3.50 DΔ					
Transitions [®] Signature GEN8 ⁽⁴⁾	60	-10.00	+6.00	-6.00	1.00 to 3.00
	65	-8.75	+6.00	-6.00	
	70	-7.75	+6.00	-4.00	
	75	-6.75	+6.00	-4.00	
Xperio [®]	60	-10.00	+8.00	-6.00	1.00 to 3.00
	65	-10.00	+8.00	-6.00	
	70	-5.00	+6.00	-6.00	
	78	-4.00	+5.00	-6.00	

Flat Top 35

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Bifocal				
	60	-8.00	+8.50	-5.00	1.00 to 3.00
	65	-8.00	+8.50	-5.00	
	70	-6.00	+8.50	-4.00	
	76	-6.00	+6.00	-3.00	
Combined 3.50 DΔ					

Round Segment 22

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Bifocal				
	67	-6.00	+18.00	-6.00	2.00 to 4.00
Combined 4.50 DΔ					

Out of range: Clear available up to +20.00 / Cyl -10.00 Combined 4.00 DΔ

TriFocal Range

Flat Top 28

1.59	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Transitions [®] Signature GEN8	Trifocal				
	82	-5.00	+7.50	-4.00	1.50 to 3.00
Combined 3.50 DΔ					
1.50 ⁽⁴⁾	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Trifocal				
	76	-8.00	+7.00	-6.00	1.50 to 3.00
Combined 3.50 DΔ					

Flat Top 35

1.50	Ø	Range	Cyl Max	ADD	Prism ⁽³⁾
Clear	Trifocal				
	75	-8.00	+7.00	-6.00	1.50 to 3.00
Combined 3.50 DΔ					

(1) Crizal[®] Easy Pro coating is not available in 1.67 index.

(2) Crizal[®] Sun XProtect[™] is only available on tinted and polarised lenses (not with mirror coatings).

(3) Prism values are indicative only and may reduce depending on powers, prism direction and effective diameter.

(4) Only available with Hard Coat and generic multicat.

Note for out of range:

Sph & Cyl values are indicative only, may reduce depending on combined powers, frame size and Pd details.

Prism values are indicative only and may reduce depending on Powers, Add, Prism direction and effective diameter.

Lens will be produced as a Biconcave in size 65 mm with an aperture 30 to 35 mm.



Essilor Exceptio on demand special range:

Unique and custom-made solutions for extreme ametropia or challenging occupational vision environments. These lenses are specially made in France, with longer lead times. Please speak to our customer service team about specific product lead times.


Range					
Single Vision					
Index	Product	Range		Cyl Max	ADD
1.50	SV Exceptio	-16.00	+6.50	-9.00	N/A
	SV Exceptio Aspheric Lenticular 1.50	-21.00	+20.00	-9.00	
	SV Exceptio Divemask Plan Concave 1.50	-10.00	-	-6.00	
	SV Exceptio Divemask Plan Convex 1.50	-	+10.00	-6.00	
	SV Exceptio Facet 1.50	-21.00	+6.50	-7.00	
	SV Exceptio Lenticular 1.50	+14.00	+34.00	-6.00	
	SV Exceptio Max Lenticular 1.50	-19.00	+20.00	-3.00	
	SV Exceptio Omega 1.50	+6.50	+20.00	-9.00	
1.67	SV Exceptio	-20.00	+20.00	-20.00	
	SV Exceptio 1.67 Transitions Grey / Brown	-20.00	+9.00	-10.00	
	SV Exceptio Facet 1.67	-21.00	+6.50	-9.00	
	SV Exceptio Max Lenticular 1.67	-30.00	+6.00	-6.00	
	SV Exceptio Max Lenticular 1.67 Transitions Grey / Brown	-20.00	+6.00	-6.00	
Bifocal					
Index	Product	Range		Cyl Max	ADD
1.50	BF Exceptio DD (UP/DN) 1.50	-6.50	+5.75	-6.00	0.75 - 3.50
	BF Exceptio RS Amblopic Hiad 1.50	-16.00	+4.00	-6.00	0.75 - 32.00
	BF Exceptio RS 22 Aspheric Lenticular 1.50	6.00	+18.00	-6.00	2.00 - 4.00
	BF Exceptio RS 38 1.50	-16.00	+15.00	-9.00	0.75 - 5.00
	BF Exceptio RS 22 Telemil 1.50	-16.00	+20.00	-9.00	0.75 - 3.50
	BF Exceptio Telex 1.50 Transitions Grey/Brown	-12.00	+7.00	-4.00	0.75 - 3.00
Progressive Lenses					
Index	Product	Range		Cyl Max	ADD
1.50	Varilux Exceptio 1.50	-21.00	+6.50	-9.00	0.75 - 3.00
	Varilux Exceptio 1.50 Transitions Grey / Brown	-21.00	+6.50	-9.00	
	Varilux Exceptio Omega 1.50	-21.00	+6.50	-9.00	0.75 - 4.00
	Varilux Exceptio - with Segment Top 1.50	-1.00	+8.00	-4.00	0.75 - 3.50
1.67	Varilux Exceptio 1.67	-25.00	+10.00	-15.00	0.75 - 3.00
	Varilux Exceptio 1.67 Transitions Grey / Brown	-25.00	+10.00	-15.00	
	Varilux Exceptio Max Lenticular 1.67	-18.00	-6.00	-4.00	0.75 - 3.50
	Varilux Exceptio Max Lenticular 1.67 Transitions Grey / Brown	-18.00	-6.00	-4.00	

1.74



Crizal®
Sapphire™ HR

1.74 Crizal® Sapphire™ HR




1.67




Crizal®
Sapphire™ HR

1.67 Crizal® Sapphire™ HR




Crizal® + BlueUV
Sapphire™ HR Capture™

1.67 Crizal® Sapphire™ HR
with Blue UV Capture™



Crizal® + Transitions®
Sapphire™ HR Gen S™

1.67 Crizal® Sapphire™ HR
Transitions™ GEN S™ Grey




1.60




Crizal[®]
Sapphire[™] HR

1.60 Crizal[®] Sapphire[™] HR




Crizal[®] + **BlueUV**
Sapphire[™] HR + **Capture[™]**

1.60 Crizal[®] Sapphire[™] HR
with Blue UV Capture[™]



Crizal[®] + **Transitions[®]**
Sapphire[™] HR + **Gen[®] S**

1.60 Crizal[®] Sapphire[™] HR
Transitions[™] GEN S[®] Grey, Amber,
Amethyst, Emerald, Sapphire, Ruby




1.60




Crizal® Rock™ + BlueUV Capture™

1.60 Crizal® Rock™ with Blue UV Capture™



Eyezen™ Boost + Crizal® Previncia®

1.60 Eyezen™ Boost⁽¹⁾ with Crizal® Previncia®



(1) Eyezen™ available in Plus, Pro and Max

1.59






Crizal®
Sapphire™^{HR}

Crizal®
Previncia®

1.59 Hard Coat

1.59 Crizal® Sapphire™ HR

1.59 Crizal® Previncia™




1.59




Crizal® + **BlueUV**
Sapphire™ HR + **Capture™**

1.59 Crizal® Sapphire™ HR
with Blue UV Capture™




Crizal® + **BlueUV**
Rock™ + **Capture™**

1.59 Crizal® Rock™
with Blue UV Capture™



Eyezen™ + **Crizal®**
Boost + **Previncia®**


1.59 Eyezen™ Boost⁽¹⁾
with Crizal® Previncia®



(1) Eyezen™ available in Plus, Pro and Max


Transitions®
Gen®S™

1.59 Hard Coat
Transitions™ GEN S™



Crizal® + **Transitions®**
Sapphire™ HR + **Gen®S™**

1.59 Crizal® Sapphire™ HR
Transitions™ GEN S™ Grey



1.50






Crizal®
Sapphire™^{HR}

1.50 Crizal® Sapphire™ HR

Crizal®
Previncia®


1.50 Crizal® Previncia™

1.50 Hard Coat




Crizal®
Rock™

1.50 Crizal® Rock™



Crizal®
Easy^{Pro}

1.50 Crizal® Easy^{Pro}




1.50




Crizal[®] Sapphire[™] HR + **BlueUV Capture[™]**

1.50 Crizal[®] Sapphire[™] HR
with Blue UV Capture[™]




Crizal[®] Rock[™] + **BlueUV Capture[™]**

1.50 Crizal[®] Rock[™]
with Blue UV Capture[™]



Crizal[®] Easy^{Pro} + **BlueUV Capture[™]**

1.50 Crizal[®] Easy^{Pro}
with Blue UV Capture[™]




1.50




Transitions
gen#S™

1.50 Hard Coat
Transitions™ GEN S® Grey




Crizal®
Sapphire™ HR + Transitions
gen#S™

1.50 Crizal® Sapphire™ HR
Transitions™ GEN S® Grey




Crizal®
Rock™ + Transitions
gen#S™

1.50 Crizal® Rock
Transitions™ GEN S® Grey




Crizal®
EasyPro + Transitions
gen#S™

1.50 Crizal® EasyPro
Transitions™ GEN S® Grey




Eyezen
Boost + Crizal®
Prevenica®

1.50 Eyezen™ Boost⁽¹⁾
with Crizal® Prevenica®



Crizal®
Sapphire™ HR + Transitions
gen#S™

1.60 Crizal® Sapphire™ HR
Transitions™ GEN S® Amber, Amethyst,
Emerald, Sapphire, Ruby




(1) Eyezen™ available in Plus, Pro and Max

Varilux® XR Base Charts

Base matching available. One base higher or lower than the base shown in the chart (e.g. 1.50 Index, +3.50 will be made on a 8.00 base, but can be made on a 6.25 base in order to match the other lens). This is a guideline only. Customer Service can advise the recommended base curve for a given prescription and material index.
NOTE: where there are two base curves indicated for a colour, the higher value is applied to higher additions.
Customer service can advise of the exact base curve selected for your order.


1.50





A base curve is assigned for a certain prescription range on the basis of being as close as possible to “best form” for those powers. Matching base curve alters the assigned base so, in some instances, this could impact on lens performance.

1.59

1.59 XR Clear/Transitions™ Base Ranges




1.59 XR Polar Base Range





A base curve is assigned for a certain prescription range on the basis of being as close as possible to “best form” for those powers. Matching base curve alters the assigned base so, in some instances, this could impact on lens performance.

1.60

1.60 XR Clear/Transitions[™] Base Ranges



1.60 XR Polar Base Range




A base curve is assigned for a certain prescription range on the basis of being as close as possible to “best form” for those powers. Matching base curve alters the assigned base so, in some instances, this could impact on lens performance.



1.67

1.74

1.67 XR Clear/Tinted/Transitions™ Base Ranges



1.74 XR Clear/Tinted/ Transitions™ Base Ranges



A base curve is assigned for a certain prescription range on the basis of being as close as possible to “best form” for those powers. Matching base curve alters the assigned base so, in some instances, this could impact on lens performance.



Varilux® XR Series						
Varilux® XR series™	Varilux® XR Pro	+ 4D	Varilux® XR 4D	+ 4D	Varilux® XR Track	/
	Varilux® XR Fit	+ F	Varilux® XR Short	+ S	Varilux® XR Design	/

Varilux® Physio® 3.0						
Varilux® Physio® 3.0	Varilux® Physio® 3.0 Eyecode	+ E	Varilux® Physio® 3.0 Fit	+ F	Varilux® Physio® 3.0 Short	+ H
	Varilux® Physio® 3.0	+ H				

Varilux® Comfort Max						
Varilux® Comfort Max	Varilux® Comfort Max Short	+ S M	Varilux® Comfort Max	+ M	Varilux® Comfort Max Eyecode	+ E

Essilor Essentials Assent						
Essentials	Essilor Essentials Assent Short	Si S	Essilor Essentials Assent	Si		

Essilor Essentials Accord						
Essentials	Essilor Essentials Accord Short	HS	Essilor Essentials Accord	H		

Varilux® Digitime®						
Varilux® Digitime®	Varilux® Digitime Near Eyecode	+ E	Varilux® Digitime Mid Eyecode	+ E	Varilux® Digitime Far Eyecode	+ E
	Varilux® Digitime Near Fit	+ F	Varilux® Digitime Mid Fit	+ F	Varilux® Digitime Far Fit	+ F
	Varilux® Digitime Near	/	Varilux® Digitime Mid	/	Varilux® Digitime Far	/

Eyezen™ Boost						
Eyezen™ Boost	Eyezen™ Boost Plus 0.40 Eyecode	+ E	Eyezen™ Boost Pro 0.60 Eyecode	+ E	Eyezen™ Boost Max 0.85 Eyecode	+ E
	Eyezen™ Boost Plus 0.40 Fit	+ F	Eyezen™ Boost Pro 0.60 Fit	+ F	Eyezen™ Boost Max 0.85 Fit	+ F
	Eyezen™ Boost Plus 0.40	/	Eyezen™ Boost Pro 0.60	/	Eyezen™ Boost Max 0.85	/

Eyezen™ Start						
Eyezen™ Start	Eyezen™ Start	/	Eyezen™ Start Fit	+ F	Eyezen™ Start Eyecode	+ E

Crizal®



The Crizal® Range

Crizal® lenses are worn by millions of people around the world.
There is a Crizal® coating for every wearer, every need and every lifestyle.

Crizal – New Generation

		Clarity	Durability			Protection
		Anti-reflection	Scratch resistant	Anti-smudge/ easy clean	Dust repellent/ anti-static	ESPF/Back surface UV cut
		Virtually reflection-free	Protect from every day wear & tear	Protect from every day cleaning	Clean less often	Protects from reflected UV rays*
Crizal® Sapphire™^{HR}	Best in Class Transparency/ Best overall anti-reflection coating	 Multi-angular Technology™	 High-Resistance Technology™	 High Surface Density Process™	✓	
Crizal® Rock™	Best combination against scratches & smudges	✓	 High-Resistance Technology™	 High Surface Density Process™	✓	
Crizal® Easy^{Pro}	An easy to clean anti-reflection coating	✓	✓	 High Surface Density Process™	✓	

* Except 1.5 index clear = ESPF 10

Crizal With Blue Light Filter Technology

	Clarity	Durability			Protection	Blue-Violet light filter
	Anti-reflection	Scratch resistant	Anti-smudge/ easy clean	Dust repellent/ anti-static	ESPF/Back surface UV cut	Selectively filter blue-violet light
	Virtually reflection-free	Protect from every day wear & tear	Protect from every day cleaning	Clean less often	Protects from reflected UV rays*	Better contrast & comfort
Crizal® Prevencia™	✓	✓	✓	✓	25+ E-SPF	Light Scan Technology

* Except 1.5 index clear = ESPF 10



**BlueUV
Capture™**
**EyeProtect
System™**

E-SPF® Index

Crizal® coatings augment protection against the invisible and often irreversible dangers of UV light.

E-SPF® is an index rating the overall UV protection of a lens, developed by Essilor International and endorsed by third-party experts. Crizal® Sapphire™ HR comes with E-SPF 35™, the highest level of E-SPF® coating for a clear Crizal® lens⁽¹⁾.

Advanced Technology For Optimal Protection

Crizal® products are available with Eye Protect System™ or Blue UV Capture technologies.


Crizal®, combined with either of these two technologies, provides everyday protection against harmful UV rays along with blue-violet light filtering.

(1) Based on sliding angle criteria and stain release test, amongst globally known lenses brands by consumers (MSW brand tracking 2018 in 17 countries).

EyeProtect System™



Essilor's Blue-violet Filter
Lens Innovation




EyeProtect System™

Essilor's most
Imperceptible* Solution
to filter UV and Blue-Violet Light

It Is Not an
Anti-Reflection Coating

*Blue-violet light: between 400 and 455nm as stated by ISO TR 20772:2018. Feature obtained via an embedded Innovation. It is not an anti-reflective coating.

BlueUV
Capture™



BlueUV Capture™

FILTERS **20%** Blue-violet light
up to 455nm

3X More filtering of blue-violet light
than a standard prescription lens

100% UV cut up to 400nm,
including 1.5 Index



Transitions®
Gen S™



SUPERPOWER YOUR GLASSES. SUPERSMOOTH YOUR LIFE.

Transitions® GEN S™ is our perfect everyday lens. It is ultra-responsive to light¹ and offers a spectacular colour palette of 8 vibrant colours. These tints are fully clear indoors and beautifully coloured outdoors, providing endless possibilities of pairing to complement any look.



Effortless vision experience



Improved long lasting
photochromic performance
vs Transitions® GEN 8™²



Fully clear indoors



Amazingly fast



Darkens in seconds³



Blocks 100% UVA & UVB rays
Filters blue-violet light
indoors and outdoors⁴

8 BEAUTIFUL COLOURS :



Grey



Graphite Green



Brown



Sapphire



Amethyst



Amber



Emerald



NEW Ruby

(1) For polycarbonate & CR39 lenses across colours achieving 18% transmission @ 23°C in less than 30 seconds. (2) For grey CR39 & polycarbonate lenses, compared to the previous generation. (3) For polycarbonate & CR39 lenses across colours achieving 18% transmission at 23°C. (4) For polycarbonate and CR39 lenses across colours. Blue-violet light is between 400 and 455nm (ISO TR 20772:2018).

Transitions and the Transitions logo are registered trademarks of Transitions Optical Inc. used under license by Transitions Optical Limited. GEN S is a trademark of Transitions Optical Limited. ©2024 Transitions Optical Limited. Photochromic performance may vary across colours and lens material and is influenced by temperature and UV exposure.

XTRACTIVE® POLARIZED™

DYNAMIC POLARISED
PHOTOCHROMIC LENS



Frames by TALA - Lenses Transitions® Grey

XTRACTIVE® POLARIZED™

CLEAR INDOORS WITH A
HINT OF PROTECTIVE TINT



INDOOR CLARITY

UP TO 2X FASTER
FADE-BACK¹



RESPONSIVENESS

BLOCKS 100%
UVA & UVB



UV PROTECTION

DARKEST EVEN IN
HOT TEMPERATURES²



DARKNESS

ACTIVATES IN THE CAR
UP TO CATEGORY 2³



IN THE CAR

UP TO 90% POLARISATION
EFFICIENCY⁴



POLARISATION EFFICIENCY

BEST BLUE-VIOLET LIGHT
FILTER INDOORS⁵




BLUE-VIOLET
LIGHT FILTERING

- (1) Compared to Transitions Vantage. Based on tests across materials on grey lenses @ 23°C, using ISO 12312-1 standard tests across materials on grey lenses @ 23°C.
- (2) Based on tests on polycarbonate grey lenses; up to 10% darker than Transitions Vantage @ 23°C and up to 5% darker @ 35°C.
- (3) Based on tests across materials on grey lenses, achieving transmission below 45% @ 23°C behind a standard windscreen. The lens achieves a polarisation efficiency of up to 30% behind the windscreen.
- (4) Based on tests across materials on grey lenses @ 23°C, using ISO 12312-1 standard.
- (5) Filters up to 45% of blue-violet light indoors. Test carried out on grey CR39 and polycarbonate lenses with a premium anti-reflective coating. Blue-violet light is between 400 and 455nm (ISO TR 20772:2018)

XTRACTIVE® NEW GENERATION

XTRA DARKNESS WITH
XTRA LIGHT PROTECTION

DEFY
THE
BRIGHT



Frames by TALA - Lenses Transitions® Grey

Transitions™
Light
Intelligent
Lenses

XTRACTIVE® NEW GENERATION

CLEAR INDOORS WITH A
HINT OF PROTECTIVE TINT



INDOOR CLARITY

BLOCKS 100%
UVA & UVB



UV PROTECTION

BEST BLUE-VIOLET LIGHT
FILTER INDOORS⁶



BLUE-VIOLET
LIGHT FILTERING

DARKEST EVEN IN
HOT TEMPERATURES⁷



DARKNESS

DARKEST IN THE CAR⁷



IN THE CAR

UP TO 35% FASTER
FADE-BACK⁸



RESPONSIVENESS

- (6) The darkest in hot temperatures, in the car and offering the best overall blue violet light filtration across light situations* among clear to extra dark photochromic lenses. *Filtering blue-violet light (between 400 and 455nm ISO TR 20772:2018) among polycarbonate and CR39 grey lenses with a premium anti-reflective coating.: filtering (i) up to 45% indoors at 23°C, (ii) up to 64% behind the windshield (iii) up to 86% outdoors at 23°C and (iv) up to 83% outdoors at 35°C.
- (7) Clear to extra dark photochromic category. Tests across polycarbonate and 1.5 grey lenses using Transitions Optical's standard testing method. Outdoors at 35°C achieves <18%T. Behind the windscreen at 23°C achieves between 18%T and 43%T.
- (8) Compared to the previous Transitions XTRActive generation across materials tested on grey lenses fading back to 70% transmission at 23°C.

Transitions™




Transitions® XTRActive® style mirrors colour matching guide

Use this guide to find which Transitions® XTRActive® style mirrors colour is best with your brown or grey tinted lenses:

Bronze



Gold Pink




Best paired with Transitions® XTRActive New Generation® Brown Lenses

Silver




Forest Green



Violet



Ocean Blue



Best paired with Transitions® XTRActive New Generation® Grey Lenses


*Note: Transitions® XTRActive® New Generation Lenses with Mirror coatings are not suitable for night driving.

Essilor Sun Solutions



Ultimate Style & Protection Without Any Sacrifice

The Essilor® Sun Lens Range gives your patients the comfort and protection they need from harmful UV rays without sacrificing the style they want.



Xperio® sun lenses eliminate blinding reflected glare, offering greater visual comfort, excellent clarity of vision and improved performance for your patients.



With a choice of progressive or single vision lens options, SolarFlair lets your patients choose any frame, whether it's sport or fashion.



Transitions® Style Mirrors offer your patients a fashionable mirror effect with lenses that adapt to the surrounding light. They're the essence of light under control!


Xperio® Sun Lenses Can Make the Difference.

Whether it's prescription sunglasses for sport, driving, as a fashion accessory or because of sensitivity to light, XPERIO® sun lenses are available with your favourite Varilux®, Single Vision and Eyezen® lenses.

Xperio® Polarised Lenses

- Ideal for driving, outdoor activities, water sports and more.
- Improved visual acuity and reaction time in driving conditions*
- Better colour perception with contrast sensitivity improved by 75%**

Available in:



Grey


Brown

Graphite Green

Xperio® Tinted Lenses

- Reduces brightness and glare
- Improves contrast and clarity
- Provides visual comfort in the sun

Leading Tint Colours:



Grey


Brown

Green



Grey Blue

Polarised & Tint Gradients



Grey

Brown

Green

Violet

Blue

Purple

Xperio® Mirror Lenses:



Gold Pink

Silver

Bronze

Silver Shadow

Violet

Turquoise Blue

Forest Green

Ocean Blue

*Survey contrast sensitivity and reaction times with polarised and tinted lenses in simulated driving daylight conditions - Institute for Vision Research, Manhattan Vision Associates program #5331 - Association for Research in Vision and Ophthalmology 2009.

**Survey conducted on 27 wearers, comparing visual performance with polarised and non-polarised sunglasses under reflected glare by measuring contrast sensitivity - Saint-Maur R&D.



Xperio® Tint Colour Range

		Brown	Green	Grey	Grey Blue
Full Tints	85%				
	65%				
	35%				
	18%				
Gradient Tints	85%				
Fashion Tints		Orange	Yellow	Sky Blue	Sky Blue Gradient
		Satin Rose	Satin Rose Gradient		

Limited range availability, detailed in the Suntint Availability section.

Essilor Sun Solutions



Product		Index	XPERIO® Polarised		Custom Tint	XPERIO® Tints	
			Full Polarised	Polarised Gradient		Full Tint	Gradient Tint
Progressive	Varilux® XR Series™	1.74					
		1.67			⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60				⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●			⦿⦿⦿⦿	
		1.50	●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Varilux® Physio® Extensee™	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●	●●●●●●●●		⦿⦿⦿⦿	
		1.50	●●●●	●●●●●●●●	⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Varilux® Comfort Max	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●	●●●●●●●●		⦿⦿⦿⦿	
		1.50	●●●●	●●●●●●●●	⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Essilor Essentials™ Assent	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●			⦿⦿⦿⦿	
		1.50	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Essilor Essentials™ Accord	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●			⦿⦿⦿⦿	
		1.50	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Varilux® Digitime	1.67			⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60				⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59				⦿⦿⦿⦿	
		1.50			⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
Extended Focus	Single Vision	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●	●●●●●●●●		⦿⦿⦿⦿	
		1.50	●●●●	●●●●●●●●	⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Eyezen® Boost	1.67			⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60				⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59				⦿⦿⦿⦿	
		1.50			⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	SolarFlair	1.67	●●●●			⦿⦿⦿⦿	
		1.60	●●●●			⦿⦿⦿⦿	
		1.59	●●●●			⦿⦿⦿⦿	
		1.50	●●●●			⦿⦿⦿⦿	
	Stylistic	1.67	●●●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
		1.60	●●●●			⦿⦿⦿⦿	⦿⦿⦿⦿
		1.59	●●●●	●●●●●●●●		⦿⦿⦿⦿	
		1.50	●●●●	●●●●●●●●	⦿	⦿⦿⦿⦿	⦿⦿⦿⦿
	Bi-focal	1.59	●			⦿⦿⦿⦿	
		1.50	●●		⦿	⦿⦿⦿⦿	⦿⦿⦿⦿

Legend

● Grey

● Graphite Green

● Gradient Grey

● Gradient Green

● Gradient Blue

⦿ Grey Tint

● Green Tint

● Brown

⦿ Custom Tint

● Gradient Brown

● Gradient Violet

● Gradient Purple

● Brown Tint

● Grey Blue Tint




Essilor Sun Solutions



Product		Index	Fashion Tints		Mirrors
			Full Tint	Gradient Tint	Paired Grey or Green Tint
Progressive	Varilux® XR Series™	1.74			
		1.67			
		1.60			
		1.59			
		1.50			
	Varilux® Physio® Extensee™	1.67			
		1.60			
		1.59			
		1.50			
	Varilux® Comfort Max	1.67			
		1.60			
		1.59			
		1.50			
	Essilor Essentials™ Assent	1.67			
		1.60			
		1.59			
		1.50			
	Essilor Essentials™ Accord	1.67			
		1.60			
		1.59			
		1.50			
	Varilux® Digitime	1.67			
		1.60			
		1.59			
		1.50			
Extended Focus	Single Vision	1.67			
		1.60			
		1.59			
		1.50			
	Eyezen® Boost	1.67			
		1.60			
		1.59			
		1.50			
Other	SolarFlair	1.67			
		1.60			
		1.59			
		1.50			
	Stylistic	1.67			
		1.60			
		1.59			
		1.50			
	Bi-focal	1.59			
		1.50			

Legend

- Satin Rose
- Sky Blue
- Orange
- Yellow
- Gold Pink Mirror
- Silver Mirror
- Bronze Mirror
- Violet Mirror
- Turquoise Blue Mirror
- Forest Green Mirror
- Silver Shadow Mirror
- Ocean Blue Mirror



essilor

[essilor.com.au](https://www.essilor.com.au)

Essilor New South Wales Lab

Ph: 02 9714 4555

Fax: 02 9748 7552

E: nsw.cs@essilor.com.au

Essilor Victoria Lab

Ph: 03 9455 7300

Fax: 03 9455 7399

E: vic.cs@essilor.com.au

Essilor Western Australia Lab

Ph: 08 9473 5000

Fax: 08 9370 4048

E: wa.cs@essilor.com.au

Essilor Queensland Lab

Ph: 07 3866 5950

Fax: 07 3866 5960

E: qld.cs@essilor.com.au

Essilor South Australia Lab

Ph: 08 8292 4200

Fax: 08 8351 5479

E: sa.cs@essilor.com.au

Essilor Tasmania Lab

Ph: 03 6234 5044

Fax: 03 6234 8101

E: tas.cs@essilor.com.au